


## Immersion Bath Unit 48965



For steady-going and economical moistening of traces with ninhydrine, DFO or other liquid developers.

The motor speed is step less adjustable and the direction of rotation is changeable, so the jamming of thin material on the rolls is removable. Also the pressure of the rolls together is adjustable, so thick material is also processible.

The remaining developer can be drained off for storage by a flexible tube

## Technical data

<b>External dimensions:</b>		
Width	mm	550
Height	mm	335
Depth	mm	410
<b>Dimensions of the rolls:</b>		
Length	mm	320
Diameter over all	mm	40
Coating (chem. Resistant)	mm	20
Motor: with frequency converter 230 V / 180 W		
Fuse		6,3 A
Connected load		230 V, 180 W, 50 Hz
Weight		25 kg

**Article number 48965**